

HD Black Line Mirrors

Proven to Reduce Glare up to 80%¹



Instrumentation and
Procedural Products

See your patients—not your instruments. Glare interferes with visibility. Reducing this disruption to your visibility allows you to see more detail and offer more accurate patient care. HD Black Line Mirrors are engineered to optimize clinical outcomes by delivering superior visibility throughout any dental procedure. The reduced glare properties of the DLC coating and clarity of the HD mirror glass facilitates quicker and more accurate visibility intraorally.



Reduces
Glare Up
To 80%¹

BEFORE

AFTER


Diamond Like Carbon (DLC) Coating

The **durable black matte coating** reduces the glare compared to a standard metal mirror head and handle. This helps to **reduce strain and fatigue** as the user does not need to adapt their viewing position due to unwanted shine.



I wear loupes with a LED light and noticed a significant difference in the amount of glare reflecting back into my field of vision during patient care, allowing me to see more detail. I realized, at the end of my clinical day, I had less eye strain when I used the HD Black Line Mirrors.

Tami Wanless, RDH, MED



Superior Brilliance & Color

HD Mirror provides **superior brilliance and color** for quicker and more accurate visibility of the mouth.²

- 113%² reflection factor for exceptional image clarity
- 38.5%² brighter than rhodium coated mirror glass
- 50%² brighter than other front surface mirror glass

Enhanced Contrast & Visual Acuity

The **black color** of the coating creates a distinct contrast between the instrument, the tooth and/or the surrounding tissue allowing for easy identification intraorally.¹

Ergonomic Handle


The large diameter handle increases control, provides maximum comfort and reduces hand fatigue, creating an ergonomically friendly handle option.

Enhance Your Set-Ups with HD Black Line Mirrors



HD Black Line Mirrors are not just for hygiene exam procedures. Complete your kit no matter the setup!³


General Exam



Basic Hygiene



Exam Perio



Restorative



Endodontic



Surgical



PARTCODES & DESCRIPTIONS

Cone Socket Thread



MIR5DSHDB
#5 Double-sided HD Mirror, Cone Socket, Black Line, Single Pack



MIR5HDB
#5 Single-Sided HD Mirror, Cone Socket, Black Line, Single Pack

European Thread



MIR5DSHDBE
#5 Double-sided HD Mirror, European, Black Line, Single Pack



MIR5HDBE
#5 Single-Sided HD Mirror, European, Black Line, Single Pack

MH6B

Cone Socket Mirror Handle, Satin Steel, Black Line

MHE6B

European Mirror Handle, Satin Steel, Black Line



Visit us online at [HuFriedyGroup.com/Mirrors](https://www.HuFriedyGroup.com/Mirrors)

1) When comparing the Hu-Friedy DLC coated mirror head and handle to the Hu-Friedy non-coated stainless-steel mirror head and handle.
2) Data on file and available upon request. 3) HD Black Line Mirrors sold separately.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-285 /0922