

An IMS® Product

**COLOR CODE TABS**

Steam sterilize for at least  
4 minutes at 270°F/132°C or  
30 minutes at 250°F/121°C.  
Do not exceed 350°F/177°C.  
Recommend 30-minute dry  
time after sterilization cycle.

Manufactured for:  
**Hu-Friedy Mfg. Co., LLC**  
3232 N. Rockwell St.  
Chicago, IL 60618 | USA  
Hu-Friedy.com

IMS-1295N | 5 TABS/YELLOW

N1220/0318



+D6718234491D